

**STAND TRIPLE PANEL (TOTAL OF 810 5M LED -405 COLOR + 405 IR)**

Output Intensity (mW/cm2) +/- 5%	Out Put Wave Lengths	Typical Therapy Dosage	Pulses per seconds	Effective Light Surface Area
		Power output		
630nm & 840nm =>50mW / cm2	Red = 630 +/- 5nm	45 Joule/cm2 (15 Minutes)	0 to 5000 pps	96.86"
530nm & 840nm =>53mW / cm2	Green = 530 +/-5nm	90 Joule/cm2 (30 Minutes)		{[5.25"x 6.15"] x 3}
445nm & 840nm =>53mW / cm2	Blue = 445 +/-5 nm	135 Joule/cm2 (45 Minutes)		
575nm & 840nm =>52mW / cm2	IR = 840 +/- 10nm			
630nm => 30mW / cm2	Amber = Δ 575 +/- 7nm			
530nm => 35mW / cm2				
445nm => 32mW / cm2				
575nm => 33mW / cm2				
840nm => 45mW / cm2				